HOOFCARE LARGE ANIMAL INDEX AT BACK. Please use product codes when ordering

CATTLE SHOOF



Cattle-Shoof

all models		\$27.50	
SIZE	COLOUR CODE	MODELS AVAILABLE	ORDER CODE
X-SMALL	PINK	Standard sole	207 422
SMALL	RED	Standard sole	207 428
MEDIUM	YELLOW	Standard sole Elevated left side Elevated right side	207 434 207 432 207 433
LARGE	BLUE	Standard sole Elevated left side Elevated right side	207 419 207 417 207 418

"Shoof" is a shoe-for-a-hoof. This product is used for the treatment of lameness in cattle.

The Shoof is purposefully constructed in a very open manner. This allows the hoof to "breathe" which is essential for a fast recovery. The medication and bandage supplied with the Shoof are also essential for successful treatment.

The Shoof copper-sulphate medication is effective in curing and hardening the hoof-horn, even in the worst conditions. The special Shoof bandage supplied acts like a filter, keeping dirt out and the copper-sulphate in, even if the animal is walking in very bad conditions (eg: dairy farm roads and yards, and muddy paddocks).

Cattle Shoof can be re-used many times with our Shoof Refill-Kit components. Non-standard bandage or medication may cause failure of Shoof. Each Cattle Shoof will last years of re-use.

Refill-Kit

complete each

207 429 \$3.50

The Cattle Shoof Refill-Kit contains the following: 1. Shoof copper-sulphate medication sachet 2. Shoof special conforming 10cm elastic bandage

3. Shoof replacement tie-string These items are essential for successful re-use of your Shoof. They are also very convenient and low-

your Shoof. They are also very convenient and lowcost. Non standard items can cause Cattle Shoof failure.



Book "Cattle Lameness and Hoofcare" 200 837 \$95.95

Written by the well known UK hoofcare specialist Roger Blowey, this is a fully illustrated guide to the causes and prevention of cattle lameness. Also contains comprehensive details of hoof anatomy and trimming techniques. Colour photos. 208 pages.







Walkease™

Starter Kit complete		
small 10-pk	212 088	\$149.00
medium 10-pk	212 089	\$149.00
large 10-pk	212 094	\$149.00
mixed sizes 10-pk	212 102	\$149.00
Blocks only		
small 10-pk	212 097	\$79.00
medium 10-pk	212 099	\$79.00
large 10-pk	212 101	\$79.00
Glue only		
		+

20ml each 211 648 \$39.00 Walkease (say'walkies') is the new hoof repair block from Shoof International Ltd. Invented and developed in-house by Mr Peter Griffiths, Shoof's product development manager. Walkease has been created as an early-intervention hoof claw prosthetic. It is not intended to replace other conventional claw shoes currently on the market. The Walkease block is made from flexible EVA material, that is soothing for the cow to walk on. The greatest advantage of the Walkease is its speed and ease of application. The Starter Kit comprises 10 Walkease blocks, a tube of special Walkease 'superglue', a Surform hoof rasp (203 624) and 10 pairs of latex gloves, all packed in a useful resealable plastic container. For further use extra 10-packs of blocks and extra glue are available. A practical 'first-purchase' would be the 'mixed sizes' starter kit (212 102). This has blocks as follows: 3×3 small, 3×3 medium and 4×3 large. This will allow the user to determine the preferred block size for their own herd.

Walkease blocks can be used either way up, so suit both left and right claws. The block should stay on from 7-14 days, which is generally long enough for the treated claw to have made a recovery. For more seriously affected claws that have had invasive treatment, we recommend use of our other listed claw and hoof shoes.

Walkease is covered by design and patent applications.







Walkease[™]

Procedure for use is:

- 1. Treat the damaged claw as required.
- 2. Use the surform rasp to clean and level the healthy claw.
- 3. Apply glue sparingly to both the hoof and the walkease block. (Max 2ml).
- 4. Press the block onto the hoof for about 10 seconds.
- 5. Immediately release the leg for the cow to put weight on the Walkease. (This part is most important)



Hoof Rasp Surform

complete 203 624 \$13.50 repl't blade 203 625 \$6.95

The Surform hoof rasp uses the same blade and cutting action as the well-known Surform wood rasps. Blade is 40mm x 60mm, mounted in a plastic handle. Blade is simply pressed out to replace. A very useful tool for hoof repair work. Length 18cm o.a.

Copyright © 2009. All rights reserved